

RECTAL CANNULA



**VECTEC®**  
SINGLE USE  
LAPAROSCOPIC  
DEVICES





## Single use cannula

- Allows excellent recto sigmoid exposure.
- Facilitates double safety testing.
  - 1- Air test
  - 2- Liquid dye test



Reference	Designation	Shaft diameter	Tip diameter	Tip length	No. unit per box
RCC 320	Canule rectale	4 mm	32 mm	50 mm	10

### Product description :

- A distal cannula with:
  - 1- an ergonomic and atraumatic tip for safe progression into the rectal and sigmoid cavities.
  - 2- a flat rear part to the tip good for safety tests (air/liquid).
- A small-diameter stainless steel shaft which is flexible to adapt to the patient's anatomy.
- An ergonomic handpiece.
- T tubing:
  - 1- a length of small diameter silicone tubing with a stop-cock valve and luerlock for air/liquid injection, to be located within the sterile zone.
  - 2- a length of large diameter tubing with a red clamp for waste.



MD Class IIa

See instructions for use



Manufactured by:

Péters Surgical  
 42 rue Benoit Frachon  
 93013 Bobigny Cedex - FRANCE  
 Tel: +33 (0)1 48 10 62 62  
 Fax: +33 (0)1 48 91 22 99  
[commercial@peters-surgical.com](mailto:commercial@peters-surgical.com)  
[www.peters-surgical.com](http://www.peters-surgical.com)