


UTERINE MANIPULATOR

VECTEC[®]

SINGLE USE
LAPAROSCOPIC
DEVICES





TOTAL LAPAROSCOPIC HYSTERECTOMY



- Optimal cervix presentation
- Insulated rotating tissue presentation valve
- Independent removable tightness balloon, keeping pneumoperitoneum during colpotomy and suturing
- Easy uterine mobility in all directions
- Light-weight for an easy operator handling (180g)
- Fully assembled for immediate use
- Latex & Phtalates free

Reference	Packaging Includes	Number of units per box
MAUT 60	<ul style="list-style-type: none"> - Uterine manipulator - Cervix dilator - Balloon plug 	10



MD Class IIa
Sterilized by ETO

See instructions for use

CE
0459

G
MED NF EN ISO 13485

UNION EUROPÉENNE
en AUVERGNE

Manufactured by:
Péters Surgical
42 rue Benoit Frachon
93013 Bobigny Cedex - FRANCE

Tel: +33 (0)1 48 10 62 62
Fax: +33 (0)1 48 10 22 99
commercial@peters-surgical.com
www.peters-surgical.com