


S-CURVE URETHRAL DILATORS

AVAILABLE WITH OR WITHOUT HYDROPHILIC COATING

PRODUCT DETAILS

- Designed for ease of passage and reduce trauma by supporting the natural curvature of the male urethra.
- Smooth hydrophilic surface for effortless dilator advancement across tight strictures.
- Embossed sizes allow for quick identification and size selection.
- Available in singles or sets.



S-Curve Urethral Dilators				
Product Code	Size (Fr)	Length (cm)	Hydrophilic Coating	Set/Single
SPUD0220	8-20	37	WITHOUT	SET
SPUD0420	8-20	37	WITH	SET
SPUD0426	8-26	37	WITH	SET
SPUD0122S	22	37	WITHOUT	SINGLE
SPUD0124S	24	37	WITHOUT	SINGLE
SPUD0126S	26	37	WITHOUT	SINGLE
SPUD0322S	22	37	WITH	SINGLE
SPUD0324S	24	37	WITH	SINGLE
SPUD0326S	26	37	WITH	SINGLE

