


Single Use Flexible Hysteroscope

Ingenuity and Insight



BioSpectrum Ltd

Unit 9 Halifax Court,
Fernwood Business Park,
Newark, Nottinghamshire,
NG24 3JP

Contact us for more information

<http://www.bio-spectrum.co.uk/>
 01636 705 084
customers@bio-spectrum.co.uk
https://twitter.com/Bio_Spectrum
<https://www.linkedin.com/company/bio-spectrumuk>
<https://www.youtube.com/channel/UCBtwuYGwGDZld8EA-hWLSXA>

Follow us





Single Use Flexible Hysteroscope

PRODUCT SPECIFICATION

Mould	SHV-P1	SHV-P2
Outer Diameter of insertion Section (00)	2.8mm	5.2mm
Working length	240mm	290mm
Working Channel (ID)	1.2mm	2.2mm
Field of View	120°	
Depth of Field	0.5-120mm	
Angulation Range	Up: 130, Down: 130°	

BioSpectrum Ltd

Unit 9 Halifax Court,
Fernwood Business Park,
Newark, Nottinghamshire,
NG24 3JP

Contact us for more information

<http://www.bio-spectrum.co.uk/>
 01636 705 084
customers@bio-spectrum.co.uk
https://twitter.com/Bio_Spectrum
<https://www.linkedin.com/company/bio-spectrumuk>
<https://www.youtube.com/channel/UCBtwuYGwGDZld8EA-hWLSXA>

Follow us





Single-Use Semi-Flexible Hysteroscope

Single-Use Semi-Flexible Hysteroscope

- Anatomical design of the insertion section gives the scope a more comfortable experience for patients.
- Design has independent large water inlet channel as well as an outlet channel. A large working channel, which makes the performance more efficient and the operative field clearer.
- Horizontal operation with an ergonomic design

Applications

- Hysteroscopy
- Reproductive Aids in Reproductive Medicine
- Small Uterine Fibroids
- IUD Removal

BioSpectrum Ltd

Unit 9 Halifax Court, Fernwood
Business Park, Newark,
Nottinghamshire,
NG24 3JP

Contact us for more information

<http://www.bio-spectrum.co.uk/>
 01636 705 084
customers@bio-spectrum.co.uk
https://twitter.com/Bio_Spectrum
<https://www.linkedin.com/company/bio-spectrumuk>
<https://www.youtube.com/channel/UCBtwuYGwGDZld8EA-hWLSXA>

Follow us





Single-Use Semi-Flexible Hysteroscope

SPECIFICATION	
Model	SHV-S1
Outer Diameter(OD)	4.5mm
Working Channel (ID)	2.35mm
Working length	220mm
Total Length of Working Channel	380mm
Bending Angle	Up: 90°, down: 90°
Direction of View	0°
Field of View	120°
Depth of Field	2-50mm
Endoscope illumination	LED light is placed at the distal end of the endoscope

BioSpectrum Ltd


Unit 9 Halifax Court, Fernwood
Business Park, Newark,
Nottinghamshire,
NG24 3JP

Contact us for more information

<http://www.bio-spectrum.co.uk/>
 01636 705 084
customers@bio-spectrum.co.uk
https://twitter.com/Bio_Spectrum
<https://www.linkedin.com/company/bio-spectrumuk>
<https://www.youtube.com/channel/UCBtwuYGwGDZld8EA-hWLSXA>

Follow us





Single-Use Rigid Hysteroscope

- Smaller OD reduces the number of dilatations and thus eases patients discomfort by shortening operative time.
- Synchronized rotation $\pm 170^\circ$ function improves the flexibility precision, safety and efficiency of surgery and brings a better surgical experience for doctors and patients.
- Easy to operate, ergonomic design which supports horizontal and vertical grips.

Applications

- Endometrial Polyps
- Fibroids
- Hysteroscopy
- Tissue Biopsy
- Abnormal Endometrial Thickening
- Abnormal IUD & Removal IUD
- Uterine Adhesions
- Mediastinal Uterus
- Endometrial Excision

BioSpectrum Ltd


Unit 9 Halifax Court,
Fernwood Business Park,
Newark, Nottinghamshire,
NG24 3JP

Contact us for more information

<http://www.bio-spectrum.co.uk/>
 01636 705 084
customers@bio-spectrum.co.uk
https://twitter.com/Bio_Spectrum
<https://www.linkedin.com/company/bio-spectrumuk>
<https://www.youtube.com/channel/UCBtwuYGwGDZld8EA-hWLSXA>

Follow us





Single-Use Rigid Hysteroscope

SPECIFICATION		
Model	SHV - 4800	SHV-6800
Outer Diameter(OD)	Ø4.8mm	Ø6.8mm
Working Channel (ID)	1.90mm	4.2mm
Working length	190mm	
Total Length	350mm	
Rotation Angle	170° in both directions	
Field of View	140°	
Depth of Field	3-50mm	
Endoscope illumination	LED light is placed at the distal end of the endoscope	

BioSpectrum Ltd

 Unit 9 Halifax Court,
 Fernwood Business Park,
 Newark, Nottinghamshire,
 NG24 3JP

Contact us for more information

 <http://www.bio-spectrum.co.uk/>
 01636 705 084
 customers@bio-spectrum.co.uk
 https://twitter.com/Bio_Spectrum
 <https://www.linkedin.com/company/bio-spectrumuk>
 <https://www.youtube.com/channel/UCBtwuYGwGDZld8EA-hWLSXA>

Follow us

