



BRIVIEW Disposable Choledochoscope

Creating a difference in healthcare



For More Information
www.bio-spectrum.com

01636 705 084

Transcholedochal LCBDE: CHV - 100/100T/100D-E


9.9Fr
3.3mm

16Fr
5.3mm


Inner Channel Diameter: 2.6mm
Easier handling for big stones with
2.3mm instruments

Channel 1: Instruments

Light-guiding Fibre Material:
Enhance brightness and improve
image quality



PTCS/Transcystic LCBDE: CHV-US100/100/100D-H



Disposable Choledochoscope

MODE	CHV-US110-H	CHV-US100D-H	CHV-US100D-E	CHV-US100-H	CHV-100-E	CHV-100T-E
Field of view	110°	110°	110°	110°	110°	110°
Depth of view	5-50mm	5-50mm	5-50mm	5-50mm	5-50mm	5-50mm
Deflection up/down	120°/120°	120°/120°	160°/130°	120°/120°	160°/130°	160°/130°
Insertion dia.	2.5mm	3.3mm	5.3mm	3.0mm	4.8mm	5.3mm
Working length	600mm	600mm	380mm	600mm	380mm	380mm
Working Channel	1.2mm	1.2mm/ 0.6mm	2.2mm/ 1.0mm	1.2mm	2.2mm	2.6mm
Rotation	Y	N	N	Y	Y	Y
Lockable Deflection	Y	Y	Y	Y	Y	Y

PL-3000

Model	PL-3000
Power	Voltage 100-240 VAC Frequency 50/60Hz Rated Input 56 VA
Illumination	Lamp 2W LED Lamp Life 30000 H Colour Temperature 3000K-7000K
Functions	White balance, Luminance control (+/-), Sharpness, Gain, GAMA, IRIS, IRIS area, Enhancement, Denoise, 3D Denoise, Zoom, Freeze, Capture, Record Better, USB Storage.
Output	SDI*2, DVI*1, S-Video*1, CVBS*1, PS/2*1
Storage	USB 3.0
Dimensions	L:350mm W:355mm H:110mm
Weight	4.2Kg



Adapter

Compatible with previous
generation processor PL-2100

Accessories



Model Number	Product Name	OD(MM)	Working Length (MM)
FTL-BF-10-250-O-C-R	Disposable Biopsy Forceps	1	2500
FTL-GF-10-100-15-FFA	Retrieval Basket	1	1000
FTL-GF-10-120-15-FFA	Retrieval Basket	1	1200
FTL-GF-10-250-15-FFA	Disposable Grasping Forceps	1	2500
FTL-GF-18-100-20-FAA	Disposable Grasping Forceps	1	1000
FTL-GF-18-100-20-FAB	Retrieval Basket	1.8	1000

Contact & Visit Us



www.bio-spectrum.com



01636 705 084



customers@bio-spectrum.co.uk



**Unit 9 Halifax Court
Fernwood Business Park
Newark
Nottingham
NG24 3JP**

