



BioSpectrum

BioSpectrum Genius Morcellator

Creating a difference in healthcare



For More Information
www.bio-spectrum.com

01636 705 084

Key Features

The Simai Genius Morcellator supplied by BioSpectrum Ltd has proven to be a “game changer” in terms of speed and efficiency. With key features including:

- Lightweight ergonomic handle design
- Adjustable blade speed to return between 13-20 grams per minute
- Fast consistent suction ensuring tissue retention throughout
- Single foot pedal offering 2 stage operation - suction, suction and morcellation
- High quality and durable cutting blades available in single use or reusable option
- Large capacity collection canisters up to 8 litres
- 12L Single Use Closed Fluid System also available



Key Specification

Blade Length	360mm length, 4.8mm diameter
Blade Options	Single Use or Reusable
Suction Pressure	0-700 mmGh (0-90 Kpa)
Collection Cannisters	2 x 4 Litre Collection Cannisters
Tissue Collection	Yes contained in cannister 1



With a performance of between 13-20 grams per minute **the Genius Morcellator** really stands out from the competition. Well built with a dedicated trolley, the system is designed to be compact and elegant offering trouble free performance.

HoLEP Instrument Set

01. HD Cystoscope



- Ø 2.9mm , autoclavable, 302mm WL
- 30° / 12° viewing angle
- Probe channel: 1.2mm
- Code: AA-003-038 / AA-003-044

02. Working Element



- Material Titanium
- Code: AG-002-303
- Probe channel: 1.2mm
- Passive cutting, handle closed, 19/22Fr.

03. Outer Sheath



- With perforation, rotatable, 22Fr.
- Code: AG-002-302

04. Inner Sheath



- With ceramic isolation, rotatable
- Code: AG-002-301

05. Obturator



- Standard, 19Fr
- Code: AA-003-011

06. Viewing Obturator



- 19 Fr.
- Code: AA-003-012

HoLEP Instrument Set

07. Nephroscope

- 45° angled ocular Nephroscope
- 250mm length
- Ø 7.45mm, suction channel Ø 1.5mm, working channel Ø 4.7mm
- Code: AA-006-031



08. Outer Sheath (Second Type)

- With perforation, rotatable, 22 Fr.
- Code: AG-002-302



09. Accessories

- Adapter for laser probes - AG-100-006
- Touny-Borst Adapter - AG-002-073
- Sealing Cap for Nephroscope - AA-003-044

Contact & Visit Us



www.bio-spectrum.com



01636 705 084



customers@bio-spectrum.co.uk



Unit 9 Halifax Court
Fernwood Business Park
Newark
Nottingham
NG24 3JP

