



Urofino Single Use Digital Ureteroscopes

Creating a difference in healthcare




For More Information
www.bio-spectrum.com

01636 705 084

Introducing the next generation of digital ureteroscopes from BioSpectrum...

- Super-slim 6.6Fr Suction Ureteroscope
- Slim 7.5Fr Suction Ureteroscope
- Slim 7.5Fr Standard Ureteroscope
- 10Fr Dual Lumen Suction Ureteroscope



Specification

Product Code	UV-US120-H = Super-slim 6.6Fr Urofino with suction UV-US130-H = Slim 7.5Fr Urofino with suction UV-US110-H = Slim 7.5Fr Urofino Standard Scope UV-130D-H = 10Fr Urofino Dual Lumen Scope
Image resolution	160K
Deflection of distal tip	Up 275° Down 275°
Field of View	110° Depth of View 5-50mm
Working Channel	Super-slim 6.6Fr Urofino = 3Fr Slim 7.5Fr Urofino = 3.6Fr 10Fr Urofino Dual Lumen = 3.6Fr and 1.8FR



Specification

- Super-slim 6.6Fr, slim 7.5Fr and 10Fr outer diameter with tapered tip
- 3Fr, 3.6Fr and 3.6Fr/1.8Fr (dual lumen) working channels
- Suction feature activated by button on handle
- Locking mechanism to lock scope into position
- 180 degree independent shaft only rotation
- 160k picture resolution
- Two light lenses ensure a bright and detailed image
- Standard or US deflection available with 275 degree up and down
- Tapered distal tip for atraumatic insertion
- Lightweight and ergonomic handle design
- Portable screen can be used as standalone unit or plugged directly into your stacking system.



Specification

Safety	A new sterile scope for every patient, reducing the risk of infection
Consistency	Predictable clinical performance every time - no issue with wear and tear 24/7
Availability	No need to wait for reprocessing of reusable scopes
Cost transparency	No service or leasing agreements
Maintenance Free	No degradation, no reprocessing and no repair costs
Compatibility	Compatible with most commonly used diagnostic and therapeutic tools
Dedicated image processor	160k resolution, dedicated light source and compatible with all Seegen scopes
Additional scopes available	Cystoscopes sizes - 14Fr and 16Fr ERCP Briview Cholangioscope Choledochoscope size - 3mm / 4.7mm / 5.3mm

Contact & Visit Us



www.bio-spectrum.com



01636 705 084



customers@bio-spectrum.co.uk



Unit 9 Halifax Court
Fernwood Business Park
Newark
Nottingham
NG24 3JP

