

HoLEP Instrument Set

01. HD Cystoscope



- Ø 2.9mm , autoclavable, 302mm WL
- 30° / 12° viewing angle
- Probe channel: 1.2mm
- Code: AA-003-038 / AA-003-044

02. Working Element



- Material Titanium
- Code: AG-002-303
- Probe channel: 1.2mm
- Passive cutting, handle closed, 19/22Fr.

03. Outer Sheath



- With perforation, rotatable, 22Fr.
- Code: AG-002-302

04. Inner Sheath



- With ceramic isolation, rotatable
- Code: AG-002-301

05. Obturator



- Standard, 19Fr
- Code: AA-003-011

06. Viewing Obturator



- 19 Fr.
- Code: AA-003-012

HoLEP Instrument Set

07. Nephroscope

- 45° angled ocular Nephroscope
- 250mm length
- Ø 7.45mm, suction channel Ø 1.5mm, working channel Ø 4.7mm
- Code: AA-006-031



08. Outer Sheath (Second Type)

- With perforation, rotatable, 22 Fr.
- Code: AG-002-302



09. Accessories

- Adapter for laser probes - AG-100-006
- Touny-Borst Adapter - AG-002-073
- Sealing Cap for Nephroscope - AA-003-044